



Disc Disease Solutions, Inc.

DDS MAX

CERVICAL TRACTION COLLAR



DDS MAX

CERVICAL TRACTION COLLAR

The **DDS MAX** is the superior approach for treating neck pain. The ambulatory traction collar is soft, comfortable, lightweight and easy-to-use. When inflated, the advance longitudinal expansion system of our air bladder offers amazing cranial lift, elongating and unloading the cervical spine significantly reducing disc pressure.

The process is simple: by increasing the intervertebral space, pressure applied on the nerve root is relieved, thereby, eliminating pain while assisting active-rehabilitation.

Spinal decompression reveals promising results for the future, and effective management of patients with all types of disc disease and spinal injuries.

The **DDS MAX** consists of a unique and patented, expandable Traction Collar (rear section), Fixed Chin Collar (front section) and Hand Air Pump.

1 TRACTION COLLAR

The unique and patented expansion belt design provides traction effect.

2 FIXED CHIN COLLAR

The ergonomically designed contoured frontal collar features soft cushions to provide comfort, support and stability

3 ONE-TOUCH VELCRO TYPE CONNECTORS

Provides simple, easy-to-use fixed-points for traction effect. Stabilizes pressure on the collar bone and chest area.

4 FLEXIBLE SUPPORT PANEL

Additional support provided by the flexible panel comfortably reinforces and upholds the integrity of the traction provided by the Cervical Traction Collar

5 AIR INJECTION NOZZLE

A simple connection system to the Hand Air Pump allows for very easy pressure adjustments and deflation.



INFLATING THE TRACTION COLLAR

Hand Air Pump



Disc Disease Solutions, Inc.

100 Commerce Way, Suite 5, Hackensack, NJ 07601 U.S.A.

Tel: (888) 495-7440 Fax: (201) 880-7799

E-mail : info@discdiseasesolutions.com

www.DDSbelt.com