



PDAC Approved for E0849

DDS Cervitrac™




The **DDS Cervitrac™** is our premier cervical traction brace, and is the most comprehensive cervical traction device for treating neck pain. It's the first device that offers the option to receive pneumatic traction while lying down in the supine position or upright for mobile treatment. Never be restricted to the confines of a doorway again!

As the first cervical traction device to offer these amazing benefits, we are proud to present our DDS Cervitrac for your consideration. The advanced longitudinal pneumatic expansion system, gently lifts from the occipital and unloads considerable pressure from the cervical column and pinched nerves. Increasing the intervertebral disc space aids in the active rehabilitation and regeneration of compressed/herniated discs, pinched nerves and multiple cervical based afflictions.

DDS Cervitrac SIZING CHART

SIZE	CIRCUMFERENCE
S	11.5~14.5 inch
L	14.6~17.5 inch



Vertical Air Pressure Technology

VAP™ Technology delivers high vertical air pressure to create traction in the affected area and to protect the vertebrae against gravity. Our uniquely designed air cells evenly distribute air pressure effectively increasing the intervertebral disc space.



Before Treatment

The herniated or bulging intervertebral disc or herniated nucleus pulposus applies pressures on the nerve root, causing pain.

After Treatment

DDS helps resume the appropriate position of the intervertebral disc, or nucleus pulposus, by increasing the intervertebral space, thereby, reducing the pressure between the lumbar/cervical spine.

FEATURES OF DDS PRODUCTS

- Drug free pain relief for patients with herniated and bulging discs, muscle spasms, sciatica, spondylolysis and spinal stenosis
- Facilitate remolding and healing by reducing load off the spine
- Ideal to wear prior to activities causing pain
- Helpful for pre-surgical pain problems
- Relaxed muscle spasms, muscle cramps and strain
- Reduces low back pain caused by weight or strenuous jobs requiring heavy lifting or long hours of driving
- Helps maintain proper posture
- Very lightweight and comfortable allowing a full range of motions
- Can be worn discreetly under clothing
- Can ease job related neck/back pain



To find out more information please call us at

888•495•7440

or email us at Info@DDSMed.com

www.DDSMed.com

**Drug-Free
Back Pain Relief!**



**PDAC/
MEDICARE
APPROVED**

Get back to what you love!

Choosing DDS, Inc. is a **Hole-in-One!**

Disc Disease Solutions products are designed to encourage individuals to be completely mobile during treatment, with our revolutionary patented Vertical Air Pressure Technology, allowing for active rehabilitation. Our braces provide spinal decompression, relieving pressure on the pinched nerve, resulting in pain relief and encouraging an active lifestyle.



www.DDSMed.com



PDAC Approved for L0631/L0648

DDS 500

The DDS 500 back brace is the first patented Spinal-Air Decompression Brace LSO with Anterior and Posterior Rigid Panels.

The most versatile back brace of its kind!



DDS, Inc. Flagship Models:


DDS 500 | DDS Double | DDS Cervitrac™

All Braces PDAC & Medicare Approved

The DDS 500 and DDS Double braces are the Flagship models of our spinal-air decompression line. These bracing models have been engineered to provide back pain relief through spinal decompression. These decompression LSO's conform to the contours of the individuals own body shape when inflated, and provide unrivaled support and comfort during use. Our belts can be discreetly worn under clothing, allowing you to perform daily activities and enabling those with an active lifestyle to enjoy daily activities pain-free.

Decompression is achieved during inflation of the brace, while it naturally pushes upward, anchoring under the first rib, and gently pushes downward on the pelvic girdle. This stretches and lifts the torso vertically, giving unparalleled traction, while diverting away weight-bearing forces and pressure within the lumbar spinal region. Unlike traditional compression braces that simply squeeze the torso and provide temporary relief, the DDS line of decompression braces provide active rehabilitation by promoting mobility and regeneration of the afflicted areas. Users can easily wear our braces during activity, allowing freedom of mobility and getting back to enjoying their life!

** All DDS, Inc. braces come with a manual pump (with Easy-read air pressure gauge), laundry bag, user manual and extension hose for easy inflation/deflation.*



DDS 500 & DDS DOUBLE SIZING CHART

SIZE	TORSO MEASUREMENT	PANEL SIZE
S	26-28 inch	S
M	29-32 inch	M
L	33-35 inch	M
XL	36-38 inch	M
2XL	39-41 inch	L
3XL	42-44 inch	L
4XL	45-48 inch	XL
5XL	49-51 inch	XL
6XL	52-55 inch	2XL
7XL	56-59 inch	2XL

*Measurement around the navel.

*Extension panel available to increase size.

*6XL & 7XL not available for DDS Double.



PDAC Approved for L0637/L0650

DDS DOUBLE

The DDS Double is the ideal solution for taller patients! Additional lateral support is provided by rigid panels.