





## FEATURES OF DDS PRODUCTS

- Pain relief for patients with herniated and bulging discs, muscle spasms, sciatica, spondylolysis and spinal stenosis
- Facilitate remolding and healing by reducing load off the spine
- OA knee pain relief
- Lessens the amount of pain for patients and can reduce pain medication
- Helpful for pre-surgical pain problems
- Relaxed muscle spasms, muscle cramps and strain
- Reduces pain caused by weight or strenuous jobs requiring heavy lifting or long hours of driving
- Helps maintain proper posture
- Very lightweight and comfortable allowing a full range of motions
- Can be worn discreetly under clothing




**VAP™ Technology** delivers high vertical air pressure to create traction in the affected area and to protect the vertebrae against gravity.



Our uniquely designed air cells evenly distribute air pressure effectively increasing the intervertebral disc space.

## PRINCIPLE OF OPERATION



### Before treatment

The herniated or bulging intervertebral disc or herniated nucleus pulposus applies pressures on the nerve root, causing pain.

### After treatment

DDS Neck and Back Braces help resume the appropriate position of the intervertebral disc or nucleus pulposus by increasing the intervertebral space thereby reducing the pressure between lumbar/cervical vertebrae.



## DDS G2 SIZING CHART

Size	Measurement
S	26" ~ 31"
M	32" ~ 35"
L	36" ~ 39"
XL	40" ~ 46"
2XL	47" ~ 53"
3XL	54" ~ 60"

\* Measurement around the navel  
\* No extension piece available

## KNEETRAC LITE SIZING CHART

Size	Thigh ( 6" above knee )	Knee ( Centered )	Calf ( 6" below knee )
S	15.5" ~ 18.5"	13"~14"	12"~14"
M	18.5" ~ 21"	14"~15"	14"~16"
L	21" ~ 23"	15"~17"	16"~18"



## MORE FEATURES




Distributed by:

### Disc Disease Solutions, Inc.

100 Commerce Way, Suite 5, Hackensack, NJ 07601 U.S.A.  
Tel: 888.495.7440 Fax: 201.880.7799  
E-mail: info@discdiseasesolutions.com

[www.DDSbelt.com](http://www.DDSbelt.com)



## Decompression Bracing

## DISC DISEASE SOLUTIONS

**Disc Disease Solutions** introduces a revolutionary approach in the treatment of lower back, neck pain and OA knee pain.

Our unique and patented Vertical Air Pressure Technology sets our products apart from all other ambulatory supports.

[www.DDSbelt.com](http://www.DDSbelt.com)



**Neck Pain?**


## DDS MAX

### CERVICAL TRACTION COLLAR

The **DDS MAX** is the superior approach for treating neck pain. The ambulatory traction collar is soft, comfortable, lightweight and easy-to-use. When inflated, the advance longitudinal expansion system of our air bladder offers amazing cranial lift, elongating and unloading the cervical spine significantly reducing disc pressure.

The process is simple: by increasing the intervertebral space, pressure applied on the nerve root is relieved, thereby, eliminating pain while assisting active-rehabilitation.

Spinal decompression reveals promising results for the future, and effective management of patients with all types of disc disease and spinal injuries.



Before Inflation



After Inflation



**Back Pain?**

## DDS G2


### DECOMPRESSION BELT

**DDS G2** Decompression Belt creates spinal decompression as its inflated with air. It anchors underneath the rib cage pushing upwards and against the hip pushing downwards.


Spinal decompression stretches the torso vertically and displaces stress away from the affected disc and nerve.

Pressure and pain levels, within the lower back, is significantly reduced which will assist in active-rehabilitation.

In addition, a newly designed extended front Velcro closure offers a higher size range accommodating individuals to fluctuate in size and allowing them to share w/ friends and family.




Before After



• The G2 Belt is flexible and does not restrict muscle movement. It prevents muscles of the lower back from atrophying or weakening

**VAP™**  
Technology  
Vertical Air Pressure Technology



**Knee Pain?**


## Kneetrac™ Lite

### DECOMPRESSION KNEE BRACE

**DDS KneeTrac Lite** is the perfect solution for on-the-go individuals with OA knee pain. The LTS Hinge gently pulls in opposite directions, decompresses and unloads the knee joint. This reduction of pressure supports and stabilizes the knee, eases pain and assists in active rehabilitation.

The Kneetrac Lite's sleek and featherweight design promotes OA Knee Pain Relief and can be discreetly worn under clothing allowing an active lifestyle.

- Drug free arthritis knee pain relief
- Relieves pressure on knee joint
- Supports and stabilizes the knee
- Hinge is placed closer to the patella and acts as a lever to displace body weight away from the knee



General Frame KneeTrac Lite



BEFORE



AFTER

